Using the RDML package

Konstantin A. Blagodatskikh, Michał Burdukiewicz, Stefan Rödiger 2017-04-29

Contents

1	Introduction to the RDML package 1.1 Philosphy of the RDML package		1 2	
	1.1 Philosphy of the RDML package			
2	Vendors and software packages supporting RDML			
3	Structure of the RDML package		2	
	3.1 Opening and observing RDML file		3	
	3.2 Copying RDML objects			
	3.3 Modifying RDML objects			
	3.4 AsTable() method			
	3.5 Getting fluorescence data		7	
	3.6 Setting fluorescence data			
	3.7 Merging RDML objects			
	3.8 Saving RDML object as RDML file			
	3.9 Creating custom functions		10	
	3.10 Creating RDML from table		12	
	3.11 Creating RDML from ABI 7500 v.2 software			
	3.12 Creating RDML from .csv files			
	3.13 Creating RDML from .xls or .xslx files			
	3.14 Functional style		14	
	3.15 Summary			
4	Creating RDML file from raw data		15	
	4.1 Initialization		15	
	4.2 Creating custom method for RDML class			
	4.3 RDML object creation			
	4.4 Setting additional information			
	4.5 Processing data			
	4.6 Export to file			
	4.7 Benchmark of the RDML generation			
	4.8 Dealing with other data formats			
	4.9 Conclusion			
5	Benchmark 2			
6	Validation of RDML files created by the RDML package		24	
Re	ferences		24	

1 Introduction to the RDML package

The RDML package was created to work with the Real-time PCR Data Markup Language (RDML) – a structured and universal data standard for exchanging quantitative PCR (qPCR) data (Lefever et al. 2009;

Jan M Ruijter et al. 2015). RDML belongs to the family of eXtensible Markup Languages (XML). It contains fluorescence data and information about the qPCR experiment. A description and the RDML schema and the RDML format are available at http://rdml.org.

The XML technology is commonly used in the **R** environment and there are several tools to manipulate XML files (Nolan and Lang 2014). When working with XML data, the workhorse package is **XML** (Lang and CRAN Team 2015). Therefore, investments development, implementation and maintenance is moderate.

We use the **R6** (Chang 2015), assertthat (Wickham 2013), plyr (Wickham 2011), dplyr (Wickham and Francois 2015), tidyr (Wickham 2014) and rlist (Ren 2015) packages.

1.1 Philosphy of the RDML package

RDML imports various data formats (CSV, XMLX) besides the RDML format. Provided that the raw data have a defined structure (as described in the vignette) the import should be done by a few clicks. The example below shows the import of amplification curve data, which were stored in a CSV file. The function rdmlEdit() was used in the RKWard IDE/GUI (Rödiger et al. 2012) (tested with version 0.6.9z+0.7.0+devel1, tested on Kubuntu 17.04; NOTE: problems were reported on systems where not the webkit component was used for the rendering of the rdmlEdit() GUI) for further processing.

Once imported enables rdmlEDIT() and other functions from the RDML package complex data visualization and processing in the R statistical computing environment.

1.2 Major functionalities of the RDML package

The public methods of the main R6 "RDML" class can be used to access and process the internally stored RDML data. These methods include:

- **new()** creates a new RDML object, empty or from a specified RDML file;
- **AsDendrogram()** plots the structure of an RDML object as a dendrogram;
- **AsTable()** represents the data contained in an RDML object (except fluorescence data) as a data.frame;
- **GetFData()** gets fluorescence data;
- **SetFData()** sets fluorescence data to an RDML object;
- AsXML() saves an RDML object as an RDML_{v.}1.2 file.

2 Vendors and software packages supporting RDML

The RDML format is supported by Bio-Rad (CFX 96 and CFX 384), Life Technologies (StepOne, ViiA7, QuantStudio) and Roche (LightCycler~96) thermo-cycler systems. In addition, several software packages (e.g., **Primer3Plus**, **RDML-Ninja**, **qBase+**) exist, which support RDML (Untergasser et al. 2007; Hellemans et al. 2007; J M Ruijter et al. 2009; Pabinger et al. 2014; Jan M Ruijter et al. 2015).

3 Structure of the RDML package

The structure of the **RDML** package mimics the RDML format and provides several R6 classes, which corresponds to RDML v1.2 format types. All major manipulations with RDML data can be done by a class called **RDML** through its public methods:

- \$new() creates new **RDML** object (empty or from specified RDML file)
- \$AsTable() represents data contained in RDML object (except fluorescent data) as data.frame.
- \$GetFData() gets fluorescent data.

- \$SetFData() sets fluorescent data to RDML object.
- \$AsDendrogram() represents structure of RDML object as dendrogram.

3.1 Opening and observing RDML file

In this section we will use the built-in RDML example file lc96_bACTXY.rdml. This file was obtained during the measurement of human DNA concentration by a *LightCycler 96* (Roche) and the *XY-Detect* kit (Syntol, Russia).

To open the lc96_bACTXY.rdml file we have to create a new RDML object with its class initializer - \$new() and the file name as parameter filename.

```
filename <- system.file("extdata/lc96_bACTXY.rdml", package = "RDML")
lc96 <- RDML$new(filename = filename)</pre>
```

Next we can check structure of our new object - lc96 by printing it.

```
lc96

#> dateMade: 2014-08-19T11:25:48

#> dateUpdated: 2014-08-27T12:06:21

#> id: [Roche Diagnostics]

#> experimenter: []

#> documentation: []

#> dye: [FAM, Hex, Texas Red, Cy5]

#> sample: [Sample 39, Sample 41, Sample 43, Sample 45, Sample 51, Sample 53, Sample 54, Sample 55, Sa

#> target: [30116ec1-44f6-4c9c-9c69-5d6f00226d4e, 69b0b5cd-591c-4012-a995-7a8b53861548, 7797a698-1b2d-1844]

#> thermalCyclingConditions: [2f78ed33-724e-4a29-97e9-92296eb868e1]

#> experiment: [ca1eb225-ecea-4793-9804-87bfbb45f81d]
```

As a result we can see field names and after::

- names of the R6 objects contained at this field after ~,
- contained values after:,
- names of list elements enclosed in [].

The fields names for all **RDML** package classes correspond to fields names of RDML types as described at http://rdml.org/files.php?v=1.2.

For the base class \mathbf{RDML} they are:

- dateMade
- dateUpdated
- id publisher and id to the RDML file.
- experimenter contact details of the experimenter.
- documentation these elements should be used if the same description applies to many samples, targets, or experiments.
- dye information about a dye.
- sample defined template solutions.
- target defined PCR reactions.
- thermalCyclingConditions cycling programs for PCR or to amplify cDNA.
- experiment

These fields can be divided by two parts:


3.1.1 Experiment field

Contains one or more experiments with fluorescence data. Fluorescence data are stored at the data level of an experiment. E.g., fluorescence data for reaction tube 45 and target bACT can be accessed with the following code:

```
fdata <-
  1c96$
    experiment $ ca1eb225 - ecea - 4793 - 9804 - 87bfbb45f81d $
    run$`65aeb1ec-b377-4ef6-b03f-92898d47488b`$
    react$ \ 45 \ $
    data$bACT$
    adp$fpoints #'adp' means amplification data points (qPCR)
head(fdata)
      cyc
               tmp
                       fluor
#> 1:
        1 68.0054 0.0782385
#> 2:
        2 68.0429 0.0753689
#> 3:
        3 68.0451 0.0736838
#> 4:
        4 68.0525 0.0723196
        5 68.0537 0.0717019
#> 5:
#> 6:
        6 68.0538 0.0714182
```

Structure of experiments can be visualized by plotting dendrogram.

lc96\$AsDendrogram()



In this dendrogram 1 we can see that our file consists of one experiment and one run. Four targets, each with two sample types (std - standard, unkn - unknown), are part of the experiment. There is only qPCR data -adp in this experiment. Ten reactions (tubes) for standard type (std) and six reaction for the unknown (unkn) type. The total number of reactions can be more than number of reactions on the plate because one tube can contain more than one target (e.g., multiplexing).

3.1.2 Additional information fields

All fields other than **experiment**. This additional information can be referenced in other parts of the RDML file. E.g., to access sample added to react 39 and get its quantity we can use code like this:

```
run$`65aeb1ec-b377-4ef6-b03f-92898d47488b`$
    react$`39`$
    sample$id
sample <- lc96$sample[[ref]]
sample$quantity$value
#> NULL
```

3.2 Copying RDML objects


R6 objects are environments, that's why simple copying results in creating reference to existing object. Then modifying of copy leads to modification of original object. To create *real* copy of object we have to use method \$clone(deep = TRUE) provided by **R6** class.

From example above we can see that modification of id2 led to modification of original object id1 but modification of cloned object id3 didn't.

3.3 Modifying RDML objects

To modify content of **RDML** objects we can use fields as setters. These setters provide type safe modification by input validation. In addition, setting lists of objects generates names of list elements.

```
# Create 'real' copy of object
experiment <- lc96$experiment$`ca1eb225-ecea-4793-9804-87bfbb45f81d`$clone(deep = TRUE)</pre>
# Try to set 'id' with wrong input type.
# Correct type 'idType' can be seen at error message.
tryCatch(experiment$id <- "exp1",</pre>
         error = function(e) print(e))
#> <simpleError in render("./vignettes/RDML.Rmd"): Assertion on 'id' failed: Must have class 'idType',</pre>
# Set 'id' with correct input type - 'idType'
experiment$id <- idType$new("exp1")</pre>
# Similar operations for 'run'
run <- experiment$run$`65aeb1ec-b377-4ef6-b03f-92898d47488b`$clone(deep = TRUE)
run$id <- idType$new("run1")</pre>
# Replace original elements with modified
experiment$run <- list(run)</pre>
lc96$experiment <- list(experiment)</pre>
lc96$AsDendrogram()
```



And we can see in 2 our modification with \$AsDendrogram() method.

3.4 AsTable() method

To get information about all fluorescence data in RDML file (type of added sample, used target, starting quantity etc.) as data.frame we can use \$AsTable() method. By default, it provides such information as:

- fdata.name aggregated name for current fluorescence data. Default pattern is position_sample_sample.type_target (e.g., D03_Sample 39_std_bACT). This pattern can be modified by name.pattern argument.
- exp.id experiment id (e.g., exp1).
- run.id run id (e.g., run1).
- react.id react (tube) id (e.g., 39).
- position react (tube) position (e.g., D03).
- sample name of the added sample (e.g., Sample 39).
- target detection target (e.g., bACT).
- target.dyeId detection dye (e.g., FAM).
- sample, type type of the added sample (e.g., std).
- adp TRUE if contains qPCR data.
- mdp TRUE if contains melting data.

To add custom columns for output *data.frame* we should pass it as named method argument with generating expression. Values of default columns can be used at custom name pattern and new columns referring to their names. Next example shows how to use \$AsTable() method with a custom name pattern and additional column.

```
else value
 }
)
# Remove row names for compact printing
rownames(tab) <- NULL
head(tab)
#>
                         fdata.name exp.id run.id react.id position
                                                                       sample
#> 1:
          D03~Sample 39~std~IPC~Cy5
                                                        39
                                    exp1
                                             run1
                                                                DO3 Sample 39
#> 2:
           DO3~Sample 39~std~X~Hex
                                                        39
                                                                DO3 Sample 39
                                     exp1
                                             run1
#> 3: D03~Sample 39~std~Y~Texas Red
                                                        39
                                      exp1
                                             run1
                                                                DO3 Sample 39
#> 4:
       DO3~Sample 39~std~bACT~FAM exp1
                                             run1
                                                        39
                                                                DO3 Sample 39
#> 5:
         DO4~Sample 39~std~IPC~Cy5
                                      exp1
                                             run1
                                                        40
                                                                DO4 Sample 39
#> 6:
           DO4~Sample 39~std~X~Hex
                                                                DO4 Sample 39
                                      exp1
                                                        40
                                             run1
#>
      target target.dyeId sample.type adp
                                             mdp
#> 1:
        IPC
                      Cy5
                                  std TRUE FALSE
#> 2:
          X
                      Hex
                                 std TRUE FALSE
#> 3:
          Y
                Texas Red
                                  std TRUE FALSE
#> 4:
       bACT
                      FAM
                                  std TRUE FALSE
#> 5:
        IPC
                                  std TRUE FALSE
                      Cy5
                                  std TRUE FALSE
#> 6:
                      Hex
```

Also, the generated data.frame is used as a query in \$GetFData() and \$SetFData() methods (see further sections).


3.5 Getting fluorescence data

We can get the fluorescence data two ways:

- direct access to data as it was described at Experiment field subsection
- using special method \$GetFData()

Advantage of \$GetFData() is that it can combine fluorescence data from whole plate to one data.frame. Major argument of this function is request, which defines fluorescence data to be got. This request is output from \$AsTable() method and can be filtered with ease by the dplyr filter() function. Also limits of cycles, output data.frame format and data type (fdata.type = 'adp' for qPCR, fdata.type = 'mdp' for melting data) can by specified (see examples below).

```
library(dplyr)
#>
#> Attaching package: 'dplyr'
#> The following objects are masked from 'package:stats':
#>
#>
       filter, lag
#> The following objects are masked from 'package:base':
#>
       intersect, setdiff, setequal, union
library(ggplot2)
# Prepare request to get only 'std' type samples
filtered.tab <- filter(tab,
                       sample.type == "std")
fdata <- lc96$GetFData(filtered.tab,</pre>
                        # long table format for usage with ggplot2
```



Our curves are not background subtrackted ass visible in the plot. To do this we use the CPP() function from the chipPCR package (Rödiger, Burdukiewicz, and Schierack 2015) as described in Rödiger et al. (2015).

```
library(chipPCR)
tab <- lc96$AsTable(</pre>
  # Custom name pattern 'position~sample~sample.type~target~run.id'
  name.pattern = paste(
             react$position,
             react$sample$id,
             private$.sample[[react$sample$id]]$type$value,
             data$tar$id,
             run$id$id, # run id added to names
             sep = "~"))
# Get all fluorescence data
fdata <- as.data.frame(lc96$GetFData(tab,
                                      # We don't need long table format for CPP()
                                      long.table = FALSE))
fdata.cpp <- cbind(cyc = fdata[, 1],</pre>
                    apply(fdata[, -1], 2,
                          function(x) CPP(fdata[, 1],
                                           x)$y))
```

Now we have preprocessed data, which we will add to our object and use during next section.


3.6 Setting fluorescence data

To set fluorescence data to RDML object we can use \$SetFData() method. It takes three arguments:

- fdata fluorescence data in long.table = FALSE format
- request output from AsTable() function, which is used as path to set data;
- fdata.type fdata.type = 'adp' for qPCR, fdata.type = 'mdp' for melting data.


Next we will set preprocessed fluorescence data to the new run $- run1_cpp$. Such subelements of RDML as experiment, run, react and data that do not exist at **RDML** object create by SetFData automaticaly (read more at Creating **RDML** from table section).

Note that colnames in fdata and fdata.name in request have to be the same!



3.7 Merging RDML objects

Merging **RDML** objects can be done by MergeRDMLs() function. It takes list of **RDML** objects and returns one **RDML** object.



3.8 Saving RDML object as RDML file

To save **RDML** object as RDML file v1.2 we can use \$AsXML() method where file.name argument is name of new RDML file. Without file.name function returns XML tree.

```
lc96$AsXML("lc96.rdml")
```

You can use RDML-ninja or the RDML validator from the RMDL consortiun to validate a RDML file created by the ${\bf RDML}$ package file.

3.9 Creating custom functions

R6 classes allow add methods to existing classes. This can be done using the \$set() method. Suppose that we decided add method to preprocess all fluorescence data and calculate Cq:

```
RDML$set("public", "CalcCq",
         function() {
           library(chipPCR)
           fdata <- as.data.frame(self$GetFData(</pre>
              self$AsTable()))
           fdata <- cbind(cyc = fdata[, 1],</pre>
                            apply(fdata[, -1],
                                  2,
                                  function(x)
                                    # Data preprocessing
                                    CPP(fdata[, 1],
                                        x)$y)
                            )
           apply(fdata[, -1], 2,
                  function(x) {
                    tryCatch(
                      # Calculate Cq
                      th.cyc(fdata[, 1], x,
                              auto = TRUE)@.Data[1],
                      error = function(e) NA)
                  })
         }
)
# Create new object with our advanced class
stepone <- RDML$new(paste0(path.package("RDML"),</pre>
                             "/extdata/", "stepone_std.rdml"))
```

And then apply our new method:


```
stepone$CalcCq()
#>
           A01_NTC_RNase P_ntc_RNase P
                                                 A02_NTC_RNase P_ntc_RNase P
#>
                                                                     12.964268
                               13.801661
           A03_NTC_RNase P_ntc_RNase P
#>
                                               A04 pop1 RNase P unkn RNase P
#>
                               14.399910
                                                                     13.529380
         A05_pop1_RNase P_unkn_RNase P
#>
                                               A06_pop1_RNase P_unkn_RNase P
#>
                               14.609048
                                                                     14.160886
#>
         A07_pop2_RNase P_unkn_RNase P
                                               A08_pop2_RNase P_unkn_RNase P
#>
                              11.839277
                                                                     10.947218
#>
         B01_pop2_RNase P_unkn_RNase P B02_STD_RNase P_10000.0_std_RNase P
#>
                               13.167509
                                                                     13.415532
   B03\_STD\_RNase\ P\_10000.0\_std\_RNase\ P\ B04\_STD\_RNase\ P\_10000.0\_std\_RNase\ P
#>
#>
                              13.363262
                                                                     14.082731
#>
    BO5_STD_RNase P_5000.0_std_RNase P
                                          BO6_STD_RNase P_5000.0_std_RNase P
#>
                               13.039183
                                                                     14.012804
#>
    B07\_STD\_RNase\ P\_5000.0\_std\_RNase\ P
                                          BO8_STD_RNase P_2500.0_std_RNase P
#>
                               9.790724
                                                                     12.981337
#>
    CO1_STD_RNase P_2500.0_std_RNase P
                                          CO2_STD_RNase P_2500.0_std_RNase P
#>
                               15.666928
                                                                     13.132126
    CO3_STD_RNase P_1250.0_std_RNase P
#>
                                          CO4_STD_RNase P_1250.0_std_RNase P
#>
                               13.905808
                                                                     14.974769
#>
    CO5_STD_RNase P_1250.0_std_RNase P
                                           CO6_STD_RNase P_625.0_std_RNase P
                              13.880906
                                                                     12.850727
```


```
#> C07\_STD\_RNase\ P\_625.0\_std\_RNase\ P C08\_STD\_RNase\ P\_625.0\_std\_RNase\ P #> 16.915642 13.706484
```

3.10 Creating RDML from table

RDML objects can be generated not only from files but from user data contained in data.frames. To do this you have to create empty **RDML** object, create data.frame, which describes data and set data by \$SetFData() method. Minimal needed information (samples, targets, dyes) will be created from data description.


```
### Create simulated data with AmpSim() from chipPCR package
# Cq for data to be generated
Cqs <- c(15, 17, 19, 21)
# PCR si,ulation will be 35 cycles
fdata <- data.frame(cyc = 1:35)</pre>
for(Cq in Cqs) {
  fdata <- cbind(fdata,
                 AmpSim(cyc = 1:35, Cq = Cq)[, 2])
}
# Set names for fluorescence curves
colnames(fdata)[2:5] <- c("c1", "c2", "c3", "c4")</pre>
# Create minimal description
descr <- data.frame(</pre>
  fdata.name = c("c1", "c2", "c3", "c4"),
  exp.id = c("exp1", "exp1", "exp1", "exp1"),
  run.id = c("run1", "run1", "run1", "run1"),
  react.id = c(1, 1, 2, 2),
  sample = c("s1", "s1", "s2", "s2"),
  sample.type = c("unkn", "unkn", "unkn", "unkn"),
  target = c("gene1", "gene2", "gene1", "gene2"),
  target.dyeId = c("FAM", "ROX", "FAM", "ROX"),
  stringsAsFactors = FALSE
)
# Create empty RDML object
sim <- RDML$new()</pre>
# Add fluorescence data
sim$SetFData(fdata, descr)
# Observe object
sim$AsDendrogram()
                                                       adp
                            gene1
                                                                2
                                          unkn
    exp1
                run1
                            gene2
                                          unkn
```






3.11 Creating RDML from ABI 7500 v.2 software

RDML objects can be created by .eds file generated by ABI 7500 v.2 or StepOne software (any of analysis inside device software have to be done!).






3.12 Creating RDML from .csv files

RDML objects can be generated by .csv file. This file have to contain first column **cyc** for qPCR data or **tmp** for melting data. Other columns – fluorescence signal.


```
filename <- system.file("extdata/from_tables/fdata.csv", package = "RDML")
csv <- RDML$new(filename = filename)
csv$AsDendrogram()</pre>
```



3.13 Creating RDML from .xls or .xslx files

RDML objects can be generated by .xls or .xslx file. This file have to contain description data.frame at sheet description and qPCR data at sheet adp and/or qPCR data at sheet mdp.

```
filename <- system.file("extdata/from_tables/table.xlsx", package = "RDML")
xslx <- RDML$new(filename = filename)
#> Warning in value[[3L]](cond): NAs introduced by coercion
#> Warning in value[[3L]](cond): NAs introduced by coercion
xslx$AsDendrogram()
```



3.14 Functional style

To provide functional programming style, which is more convenient in \mathbf{R} , the \mathbf{RDML} class methods have function wrappers:

```
obj$AsTable(...) - AsTable(obj, ...)
obj$SetFData(...) - SetFData(obj, ...)
obj$GetFData(...) - GetFData(obj, ...)
obj$AsDendrogram(...) - AsDendrogram(obj, ...)
```

3.15 Summary

RDML package provides classes and methods to work with RDML data generated by real-time quantitative PCR devices or to create RDML files from user generated data. Because classes of the **RDML** package are build with **R6** they can be modified by adding custom methods and suggest type safe usage by input validation.

4 Creating RDML file from raw data

This vignette describes important step how to create **RDML** object from user provided data and calculate Cq values with custom **RDML** class method. We use a real-world example to demonstrate the use of the **RDML** package.

For this purpose we used data from the VideoScan platform. The VideoScan platform is a highly versatile fluorescence microscope based imaging platform for the real-time analysis of biomolecular interactions. VideoScan has been extended by a heating/cooling-unit, which can be used for DNA melting curve analysis (Rödiger, Böhm, and Schimke 2013), real-time quantitative PCR (Rödiger et al. 2013) and quantitative isothermal amplification (Spiess et al. 2015). In this example we focus on qPCR data obtained from the amplification of the human gene mlc-2v (Sheikh, Lyon, and Chen 2015). The experimental details (e.g., primer sequences, amplification conditions) are described in (Rödiger et al. 2013; George et al. 2016).

In detail, this dataset contains three amplification curves collected by three runs monitored with TaqMan qPCR chemistry. Similar to other experimental platforms has the VideoScan systems an output as csv file (comma separated value). In general, csv is a format that is sufficient for an in house usage. However, the exchange with other scientists is demanding, since further information should be part of the data. These are for example information about qPCR conditions, primers, experimentalist and so on as described in the MIQE and RDML standard documentations.

All csv raw data were transferred to the **R** environment and stored as rda files in the **C54** dataset after completion of the experiment as described elsewhere (R Development Core Team 2012). The **R** environment provides numerous packages for the import of standard file formats, including formats from spread sheet applications (e.g., XLS, XLSX, ODS) (Schutten 2014, Warnes et al. (2015), Dragulescu (2014)), an ODBC connection (Ripley and 2002) 2015) or from the clipboard (read.table("clipboard")). There are several graphical user interfaces (GUIs) for **R**, which facilitate an easy import of raw data (Valero-Mora and Ledesma 2012). For example **RKWard** (Rödiger et al. 2012) has import functionality for csv, txt, spss and stata files.

Therefore, we will use the **C54** dataset directly from the chipPCR package (Rödiger, Burdukiewicz, and Schierack 2015). The amplification curves have different number of cycles and one of them contains a missing value due to a sensor dropout. Thus, this data should be preprocessed before Cq calculation.

 ${f R}$ is a language with dynamic typing. This is helpful while scripting, but can lead to problematic debugging of more complicated workflows. ${m R6}$ classes provide type-safe interfaces to set data without access to the inner structure of objects so that all imported data can be validated. This option is highly useful when creating packages at an intermediate level for other packages (e.g., such an approach does not permit set type ${\it character}$ in place of type ${\it integer}$). Furthermore, the inheritance of ${\it R6}$ objects unifies the structure of the package and streamlines the extending of its capabilities (since the whole package is written around a single base class).

4.1 Initialization

First of all we have to load all necessary packages and fluorescence data into our ${\bf R}$ session.

```
library(RDML)
library(chipPCR)
library(MBmca)
library(ggplot2)
data(C54)
```

We can see that our data are arranged in data.frame with cycle numbers in the first column and fluorescence values in the other columns. Such data structure is ready to use with the **RDML** package and does not need any conversion.

```
str(C54)
                    56 obs. of 4 variables:
#> 'data.frame':
                 0 1 2 3 4 5 6 7 8 9 ...
    $ Cycle: num
           : num
                  0.141 0.146 0.152 0.151 0.153 ...
                  0.148 0.148 0.148 0.149 0.149 ...
    $ D2
           : num
                  0.145 0.147 0.148 0.148 0.149 ...
           : num
dat <- data.frame(Cycle = rep(C54[, 1], ncol(C54) - 1),</pre>
                  Experiment = unlist(lapply(colnames(C54)[-1], function(i) rep(i, nrow(C54)))),
                  refMFi = unlist(C54[, -1]))
levels(dat[["Experiment"]]) <- c("D1 - Stock cDNA", "D2 - 1/10 cDNA", "D3 - 1/100 cDNA")</pre>
ggplot(dat, aes(x = Cycle, y = refMFi, colour = Experiment)) + geom_point()
  0.6 -
                                                                 Experiment
                                                                     D1 - Stock cDNA
                                                                     D2 - 1/10 cDNA
                                                                     D3 - 1/100 cDNA
  0.2 -
                          20
                                            40
                               Cycle
```

4.2 Creating custom method for RDML class

The preprocessing and Cq calculation steps are held within a **RDML** object. To do this we shall add custom method to the **RDML** class and use this advanced class as **RDML** object initializator (data import).

As described previously (Rödiger et al. 2015), we first preprocess the raw data with the \mathbf{CPP} function from the chipPCR package. Next, these preprocessed data are stored in new runs within the same experiment. For Cq calculation we used the $\mathbf{diffQ2}$ function from the MBmca (Rödiger, Böhm, and Schimke 2013) package. Note, there are other functions available within \mathbf{R} to determine the Cq values. The qpcR package is one of the most popular packages for such tasks.

Additional our method will have two arguments:

- vector last.cycle cycle limits for every amplification curve;
- list bg.range background range for every amplification curve.

```
RDML$set("public", "ProcessVideoScan",
         function(last.cycle,
                  bg.range) {
           # Get curves description
           tab <- as.data.frame(self$AsTable())</pre>
           # Get fluorescent point
           dat <- as.data.frame(self$GetFData(tab))</pre>
           # Give new names to runs
           # Preprocess amplification curve raw data as described in
           # Rö diger et al. (2015) Bioinformatics. Note that the dataset
           # has different length of cycle numbers and missing values (NAs).
           # The CPP function from the chipPCR package prepares the data for
           # further analysis.
           tab[, "run.id"] <- pasteO(tab[, "run.id"], "_CPP")</pre>
           # Curves will be used one by one
           for(i in 2:length(dat)) {
             # Preprocess data
             preprocessed <- CPP(dat[1:last.cycle[i - 1], 1],</pre>
                                  dat[1:last.cycle[i - 1], i],
                                  trans = TRUE,
                                  bg.range = bg.range[[i - 1]])[["y.norm"]]
             # Create data.frame with cycle number and preprocessed curve Then
             # give name to fluorescence points columns as before preprocessing
             dat CPP <- cbind(dat[1:last.cycle[i - 1], 1],</pre>
                               preprocessed)
             colnames(dat_CPP)[2] <- tab$fdata.name[i - 1]</pre>
             # Set preprocessed data with new description (new run names)
             self$SetFData(dat_CPP, tab)
             # Calculate and set Cq
             # Set Cq from second derivative maximum method as described in
             # Rö diger et al. (2015) Bioinformatics for preprocessed data.
             # The diffQ2 function from the MBmca package
             # (Reouml; diger et al. (2013), The R Journal) was used to calculate the
             # Cq values of each amplification curve.
             cq <- diffQ2(dat_CPP, inder = TRUE)[["xTm1.2.D2"]][1]</pre>
             self$experiment[[tab[i - 1, "exp.id"]]]$
               run[[tab[i - 1, "run.id"]]]$
               react[[tab[i - 1, "react.id"]]]$
               data[[tab[i - 1, "target"]]]$cq <- cq</pre>
         }, overwrite = TRUE
```

4.3 RDML object creation

Now we can create an empty **RDML** object with our advanced class and set the fluorescence data into it. For setting the data we have to create a description (metadata) of the amplification curves, which contain:

- fdata.name curve name. Have to be as columns names in data.frame with fluorescence data.
- exp.id experiment name.
- run.id run name. In our case we have three runs.
- react.id react ID. E.g., number of tube within a reaction plate.

```
sample – sample name.
type – sample type.
target – target name. E.g., gene of interest name.
target.dyeID – fluorescent dye used for detection.
```

```
# Create a data frame of metadata
description <- data.frame(</pre>
  fdata.name = c("D1", "D2", "D3"),
  exp.id = c("exp1", "exp1", "exp1"),
  run.id = c("run1", "run2", "run3"),
  react.id = c(1, 1, 1),
  sample = c("D1 - Stock cDNA", "D2 - 1/10 cDNA", "D3 - 1/100 cDNA"),
  sample.type = c("unkn", "unkn", "unkn"),
  target = c("MLC-2v", "MLC-2v", "MLC-2v"),
  target.dyeId = c("Cy5", "Cy5", "Cy5"),
  stringsAsFactors = FALSE
# Create an empty RDML object
video.scan <- RDML$new()</pre>
# Add fluorescence data and metadata to the RDML object from a given source.
# For the sake of simplicity we use the C54 dataset from the chipPCR package.
video.scan$SetFData(C54, description)
```

4.4 Setting additional information

As optional step, we can add some information about our experiment. Such as experimenter contacts, references to documentation, description of samples and targets, thermal cycling conditions, etc.

```
# Add experimentalist information
video.scan$experimenter <-</pre>
    experimenterType$new(
      idType$new("SR"),
      "Stefan",
      "Roediger",
      "stefan.roediger@b-tu.de"
    experimenterType$new(
      idType$new("CD"),
      "Claudia",
      "Deutschmann"
    )
  )
# Add a reference to documentation
video.scan$documentation <- list(</pre>
  documentationType$new(
    idType$new("Roediger et al. 2013"),
    paste("A Highly Versatile Microscope Imaging Technology Platform for the Multiplex",
          "Real-Time Detection of Biomolecules and Autoimmune Antibodies. S. Roediger,",
          "P. Schierack, A. Boehm, J. Nitschke, I. Berger, U. Froemmel, C. Schmidt, M.",
          "Ruhland, I. Schimke, D. Roggenbuck, W. Lehmann and C. Schroeder. Advances in",
```

```
"Biochemical Bioengineering/Biotechnology. 133:33-74, 2013.",
          "https://www.ncbi.nlm.nih.gov/pubmed/22437246")
  )
)
cdna <-
  cdnaSynthesisMethodType$new(
    enzyme = "SuperScript II",
    primingMethod =
      primingMethodType$new("oligo-dt"),
    dnaseTreatment = TRUE
video.scan$sample$`D1 - Stock cDNA`$description <- "Input stock cDNA was used undiluted (D1)"
video.scan$sample$`D1 - Stock cDNA`$cdnaSynthesisMethod <- cdna</pre>
video.scan$sample$`D2 - 1/10 cDNA`$description <- "1/1000 diluted in A. bidest"
video.scan$sample$`D2 - 1/10 cDNA`$cdnaSynthesisMethod <- cdna</pre>
video.scan$sample$`D3 - 1/100 cDNA`$description <- "1/1000000 diluted in A. bidest"</pre>
video.scan$sample$`D3 - 1/100 cDNA`$cdnaSynthesisMethod <- cdna</pre>
video.scan$target$`MLC-2v`$xRef <- list(</pre>
  xRefType$new("uniprot",
                "P10916")
)
video.scan$target$`MLC-2v`$sequences <-</pre>
  sequencesType$new(
    forwardPrimer <- oligoType$new(</pre>
      sequence = "ACAGGGATGGCTTCATTGAC"),
    reversePrimer <- oligoType$new(</pre>
      sequence = "ATGCGTTGAGAATGGTTTCC"),
    probe1 <- oligoType$new(</pre>
      threePrimeTag = "Atto647N",
      sequence = "CAGGGTCCGCTCCCTTAAGTTTCTCC",
      fivePrimeTag = "BHQ2")
  )
tcc <-
  thermalCyclingConditionsType$new(
    idType$new("Amplification"),
    experimenter = list(
      idReferencesType$new("SR"),
      idReferencesType$new("CD")
    ),
    step =
      list(
        stepType$new(
          nr = 1,
          temperature = temperatureType$new(95,
                                              600)
        ),
        stepType$new(
          nr = 2,
          temperature = temperatureType$new(95,
```

```
stepType$new(
          nr = 3,
          temperature = temperatureType$new(58.5,
                                             90)
        ),
        stepType$new(
          nr = 4,
          temperature = temperatureType$new(68.5,
                                             90)
        ),
        stepType$new(
          nr = 5,
          loop = loopType$new(goto = 2,
                               repeat.n = 49)
        )
      )
  )
video.scan$thermalCyclingConditions <- list(</pre>
#add description of the experiment
video.scan$experiment$exp1$description <-</pre>
  paste("The aim was to amplify MLC-2v in the VideoScan and to monitor with a",
        "hydrolysis probe for MLC-2v. The primer sequences for MLC-2v were taken",
        "from Roediger et al. (2013). The amplification was detected in solution of",
        "the 1 HCU (see Roediger et al. 2013 for details). A 20 micro L PCR reaction",
        "was composed of 250 nM primer (forward and reverse), 1x Maxima Probe qPCR",
        "Master Mix (Fermentas), 1 micro L template (MLC-2v amplification product in",
        "different dilutions), 50 nM hydrolysis probe for MLC-2v and A.",
        "bidest. During the amplification, fluorescence was measured at 59.5 degree",
        "Celsius. The Cy5 channel was used to monitor the MLC-2v specific hydrolysis",
        "probe. Input stock cDNA was used undiluted (D1). D2 was 1/1000 and D3",
        "1/1000000 diluted in A. bidest. The D1, D2, and D3 have different numbers",
        "measure points and D2 contains a missing value at cycle 37.")
video.scan$experiment$exp1$run$run1$thermalCyclingConditions <- idReferencesType$new("Amplification")</pre>
video.scan$experiment$exp1$run$run2$thermalCyclingConditions <- idReferencesType$new("Amplification")</pre>
video.scan$experiment$exp1$run$run3$thermalCyclingConditions <- idReferencesType$new("Amplification")</pre>
```


4.5 Processing data

After creating RDML object and setting data we can use our method \$ProcessVideoScan.

Then visualize the object with the \$AsDendrogram method (runs with IDs that include _CPP in names, contain


the preprocessed curves):

```
# Visualize RDML object
video.scan$AsDendrogram()
```



And plot preprocessed curves indicating the Cq values:

```
## Visualise preprocessed data with Cq values as vertical dashed lines
# Add custom column that contains the calculated Cq
tab <- video.scan$AsTable(cq = {</pre>
                 cq <- data$cq
                 if (is.null(cq) || is.na(cq))
                   NULL
                 else
                   cq
               })
# Get preprocessed data in 'long.table' format
dat <- video.scan$GetFData(tab[grep1("_CPP", tab[["run.id"]]), ],</pre>
                           long.table = TRUE)
ggplot(dat, aes(x = cyc, y = fluor)) +
  geom_line(aes(group = fdata.name, color = fdata.name),
            size = 1.2) +
  geom_vline(aes(xintercept = cq, color = fdata.name),
             linetype = "longdash") +
  scale_x_continuous("Cycle") +
  scale_y_continuous("Fluorescence") +
  ggtitle("Amplification of human MLC-2v\nVideoScan HCU")
```



However, there are further **R** packages on CRAN, Bioconductor and GitHub, which can be used with the RDML package (Rödiger et al. 2015). To deal with non-detects (McCall et al. 2014), normalization (Perkins et al. 2012), expression analysis (Dvinge and Bertone 2009, matz_no_2013) and periodicity in qPCR systems (Spiess et al. 2016).

4.6 Export to file

To save our **RDML** object as file we can use \$AsXML() method with file name as argument: video.scan\$AsXML("filename.RDML")

Note that currently due to limits of XML package it is rather slow operation!!!

A XML files can be compressed to save space. The **RMDL** package uses the **zip()** function from the **utils** package with default values to compress the RDML file. Therefore, the compression mode depends on the operating system. The **RDML** function AsXML() without an argument returns one string, which contains the XML data. This can be used to compress it with any other compression algorithm.

4.7 Benchmark of the RDML generation

The XML format becomes less efficient to work with at larger file sizes and complex structured information (e.g., matrices). Therefore, a limitation of RMDL is the process of the file generation. This characteristic is inherited by the **RDML** package. Speed is an issue in high-throughput applications. To give the user a foundation for a decision support if the implementation of the RMDL format is advisable, we performed a benchmark. We used the **microbenchmark** package (Olaf Mersmann and Friedman 2014, Radford (2014)) to assess the time required for generation of RDML files. The simulation included datasets with ranging from 1 to 1000 entries per file.

4.8 Dealing with other data formats

Researchers often collect gene expression data from more than one laboratory and and are thus required to analyze and aggregate many data sets. Adaptable data management - also known as "adaptive informatics" - is relevant in cases where data from different omics approaches and assays (e.g., flow cytometry, digital PCR, NGS) need to be merged (Baker 2012). As RDML is based on XML, it is possible to converge the files with other formats such as HDF, (Millard et al. 2011) which is supported by **R** (Fischer and Pau 2015). This enables extended data storage and analysis, but also a higher level of experimental data management.

4.9 Conclusion

In this example we used raw data from the *VideoScan* platform to demonstrate the creation of an *RDML* file via the **RDML** package. This procedure can be easily adapted to other experiments and devices and gives the researcher a tool to create a file in a standardized format for convenient data exchange. In contrast to simple *csv* files contains this file additional metadata in a defined structure. This in turn makes the data available to other scientists for reproducible research.

5 Benchmark

The RDML format may also store data from the high-throughput technologies as digital PCR or 384-well qPCR. To make the RDML package the most appropriate tool for such applications, we optimize the saving functionality of the package. Thanks to this, the RDML package is able to save in a reasonable time even large number of PCR experiments.

We benchmarked the time needed for saving RDML object with different numbers of experiments on two desktop machines:

- 1. Windows 7 x64 (build 7601) Service Pack 1 on Intel Core i5 (2,40 GHz) with 4,00 GB RAM.
- 2. Kubuntu 15.04 (Kernel 3.19.0-27-generic) Intel(R) Pentium(R) CPU 997 @ 1.60GHz with 4,00 GB RAM

In case "1." we used R version 3.2.1 (2015-06-18) and in case "2." we used R version 3.1.2 (2014-10-31).

To have the fairest comparison and take into account various factor like garbage collection, we repeated the procedure 100 times.




Table 1: Results of benchmark.

Number of experiments	Median time [s]	Operating System
1	0.0525600	Windows
8	0.1965201	Windows
20	0.4882732	Windows
30	0.7770339	Windows
32	0.8454795	Windows
50	1.4801968	Windows
96	3.8208417	Windows
100	4.0498602	Windows
300	25.6566066	Windows
384	39.9697519	Windows
765	142.1440033	Windows
1000	239.7105810	Windows
1	0.0665418	Unix
8	0.2468777	Unix
20	0.6101219	Unix
30	0.9658560	Unix
32	1.0468907	Unix
50	1.8240828	Unix
96	4.5716005	Unix
100	4.8454858	Unix
300	29.0860889	Unix
384	44.9543657	Unix
765	159.6792129	Unix
1000	266.5159046	Unix

The result above indicate that, thanks to the RDML package, even average laboratory workstation is able to save large RDML files in reasonable time (around 4 minutes for 1000 experiments).

6 Validation of RDML files created by the RDML package

The correct structure of an RDML file is an important requirement to perform reproducibel analysis. Unfortunately created files can not be properly tested now because official RDML validator rdml.org does not work from time to time. However, there is a solution to validate the files by the alternate editor - RDML-Ninja (Jan M Ruijter et al. 2015).

References

Baker, Monya. 2012. "Quantitative Data: Learning to Share." Nature Methods 9 (1): 39–41. doi:10.1038/nmeth.1815.

Chang, Winston. 2015. R6: Classes with Reference Semantics. https://CRAN.R-project.org/package=R6.

Dragulescu, Adrian A. 2014. Xlsx: Read, Write, Format Excel 2007 and Excel 97/2000/XP/2003 Files. https://CRAN.R-project.org/package=xlsx.

Dvinge, Heidi, and Paul Bertone. 2009. "HTqPCR: High-Throughput Analysis and Visualization of Quantitative Real-Time PCR Data in R." *Bioinformatics* 25 (24): 3325–6. doi:10.1093/bioinformatics/btp578.

Fischer, Bernd, and Gregoire Pau. 2015. "rhdf5: HDF5 Interface to R." R Package Version 2.10.0.

https://www.bioconductor.org/packages/release/bioc/html/rhdf5.html.

George, Sandra, Stefan Rödiger, Christian Schröder, Michael Knaut, and Jan-Heiner Küpper. 2016. "Development of Multiplex PCR Systems for Expression Profiling of Human Cardiomyocytes Induced to Proliferate by Lentivirus Transduction of Upcyte Genes." *Journal of Cellular Biotechnology* 2 (1): 35–55. doi:10.3233/JCB-15025.

Hellemans, Jan, Geert Mortier, Anne De Paepe, Frank Speleman, and Jo Vandesompele. 2007. "qBase Relative Quantification Framework and Software for Management and Automated Analysis of Real-Time Quantitative PCR Data." *Genome Biology* 8 (2): R19. doi:10.1186/gb-2007-8-2-r19.

Lang, Duncan Temple, and the CRAN Team. 2015. XML: Tools for Parsing and Generating XML Within R and S-Plus. https://CRAN.R-project.org/package=XML.

Lefever, Steve, Jan Hellemans, Filip Pattyn, Daniel R. Przybylski, Chris Taylor, René Geurts, Andreas Untergasser, and Jo Vandesompele. 2009. "RDML: Structured Language and Reporting Guidelines for Real-Time Quantitative PCR Data." *Nucleic Acids Research* 37 (7): 2065–9. https://nar.oxfordjournals.org/content/37/7/2065.

McCall, Matthew N., Helene R. McMurray, Hartmut Land, and Anthony Almudevar. 2014. "On Non-Detects in qPCR Data." *Bioinformatics* 30 (16): 2310–6. doi:10.1093/bioinformatics/btu239.

Millard, Bjorn L., Mario Niepel, Michael P. Menden, Jeremy L. Muhlich, and Peter K. Sorger. 2011. "Adaptive Informatics for Multifactorial and High-Content Biological Data." *Nature Methods* 8 (6): 487–92.

Nolan, Deborah, and Duncan Temple Lang. 2014. XML and Web Technologies for Data Sciences with R. Springer. https://www.springer.com/statistics/computational+statistics/book/978-1-4614-7899-7.

Olaf Mersmann, Rainer Hurling, Claudia Beleites, and Ari Friedman. 2014. *microbenchmark: Accurate Timing Functions*. https://CRAN.R-project.org/package=microbenchmark.

Pabinger, Stephan, Stefan Rödiger, Albert Kriegner, Klemens Vierlinger, and Andreas Weinhäusel. 2014. "A Survey of Tools for the Analysis of Quantitative PCR (qPCR) Data." *Biomolecular Detection and Quantification* 1 (1): 23–33.

Perkins, James R., John M. Dawes, Steve B. McMahon, David LH Bennett, Christine Orengo, and Matthias Kohl. 2012. "ReadqPCR and NormqPCR: R Packages for the Reading, Quality Checking and Normalisation of RT-qPCR Quantification Cycle (Cq) Data." BMC Genomics 13 (1): 296. doi:10.1186/1471-2164-13-296.

R Development Core Team. 2012. R Data Import/Export. Vienna, Austria: R Foundation for Statistical Computing. https://www.R-project.org/.

Radford, Neal. 2014. Inaccurate Results from Microbenchmark. https://web.archive.org/web/20140823221902/https://www.r-bloggers.com/inaccurate-results-from-microbenchmark/.

Ren, Kun. 2015. rlist: A Toolbox for Non-Tabular Data Manipulation. https://CRAN.R-project.org/package=rlist.

Ripley, Brian, and Michael Lapsley (1999 to Oct 2002). 2015. RODBC: ODBC Database Access. https://CRAN.R-project.org/package=RODBC.

Rödiger, Stefan, Alexander Böhm, and Ingolf Schimke. 2013. "Surface Melting Curve Analysis with R." The R Journal 5 (2): 37–53. https://journal.r-project.org/archive/2013-2/roediger-bohm-schimke.pdf.

Rödiger, Stefan, Michał Burdukiewicz, and Peter Schierack. 2015. "chipPCR: an R package to pre-process raw data of amplification curves." *Bioinformatics* 31 (17): 2900–2902. https://bioinformatics.oxfordjournals.org/content/31/17/2900.abstract.

Rödiger, Stefan, Michał Burdukiewicz, Konstantin Blagodatskikh, and Peter Schierack. 2015. "R as an Environment for the Reproducible Analysis of Dna Amplification Experiments." The R Journal 7 (2): 127–50.

Rödiger, Stefan, Thomas Friedrichsmeier, Prasenjit Kapat, and Meik Michalke. 2012. "RKWard: A Comprehensive Graphical User Interface and Integrated Development Environment for Statistical Analysis

with R." Journal of Statistical Software 49 (9): 1-34. https://www.jstatsoft.org/article/view/v049i09.

Rödiger, Stefan, Peter Schierack, Alexander Böhm, Jörg Nitschke, Ingo Berger, Ulrike Frömmel, Carsten Schmidt, et al. 2013. "A highly versatile microscope imaging technology platform for the multiplex real-time detection of biomolecules and autoimmune antibodies." *Advances in Biochemical Engineering/Biotechnology* 133: 35–74.

Ruijter, J M, C Ramakers, W M H Hoogaars, Y Karlen, O Bakker, M J B van den Hoff, and A F M Moorman. 2009. "Amplification Efficiency: Linking Baseline and Bias in the Analysis of Quantitative PCR Data." *Nucleic Acids Research* 37 (6): e45.

Ruijter, Jan M, Steve Lefever, Jasper Anckaert, Jan Hellemans, Michael W Pfaffl, Vladimir Benes, Stephen A Bustin, Jo Vandesompele, Andreas Untergasser, and others. 2015. "RDML-Ninja and RDMLdb for Standardized Exchange of qPCR Data." *BMC Bioinformatics* 16 (1). BioMed Central Ltd: 197.

Schutten, Gerrit-Jan. 2014. readODS: Read ODS Files and Puts Them into Data Frames. https://CRAN. R-project.org/package=readODS.

Sheikh, Farah, Robert C. Lyon, and Ju Chen. 2015. "Functions of Myosin Light Chain-2 (MYL2) in Cardiac Muscle and Disease." Gene 569 (1): 14–20.

Spiess, Andrej-Nikolai, Claudia Deutschmann, Michał Burdukiewicz, Ralf Himmelreich, Katharina Klat, Peter Schierack, and Stefan Rödiger. 2015. "Impact of Smoothing on Parameter Estimation in Quantitative DNA Amplification Experiments." Clinical Chemistry 61 (2): 379–88. doi:10.1373/clinchem.2014.230656.

Spiess, Andrej-Nikolai, Stefan Rödiger, Michał Burdukiewicz, Thomas Volksdorf, and Joel Tellinghuisen. 2016. "System-Specific Periodicity in Quantitative Real-Time Polymerase Chain Reaction Data Questions Threshold-Based Quantitation." Scientific Reports 6 (December): 38951. doi:10.1038/srep38951.

Untergasser, Andreas, Harm Nijveen, Xiangyu Rao, Ton Bisseling, René Geurts, and Jack A. M. Leunissen. 2007. "Primer3Plus, an Enhanced Web Interface to Primer3." *Nucleic Acids Research* 35 (suppl 2): W71–W74. doi:10.1093/nar/gkm306.

Valero-Mora, Pedro M., and Ruben Ledesma. 2012. "Graphical User Interfaces for R." *Journal of Statistical Software* 49 (1): 1–8. https://www.jstatsoft.org/article/view/v049i01.

Warnes, Gregory R., Ben Bolker, Gregor Gorjanc, Gabor Grothendieck, Ales Korosec, Thomas Lumley, Don MacQueen, Arni Magnusson, Jim Rogers, and others. 2015. gdata: Various R Programming Tools for Data Manipulation. https://CRAN.R-project.org/package=gdata.

Wickham, Hadley. 2011. "The Split-Apply-Combine Strategy for Data Analysis." *Journal of Statistical Software* 40 (1): 1–29. https://www.jstatsoft.org/article/view/v040i01.

- ——. 2013. assertthat: Easy Pre and Post Assertions. https://CRAN.R-project.org/package=assertthat.
- ——. 2014. tidyr: Easily Tidy Data with Spread() and Gather() Functions. https://CRAN.R-project.org/package=tidyr.

Wickham, Hadley, and Romain Francois. 2015. dplyr: A Grammar of Data Manipulation. https://CRAN. R-project.org/package=dplyr.